

Endosafe[®] - PTS[™] Cartridges Certificate of Analysis

Reorder Code: PTS2005F Cartridge Lot #:  8698116 Expiration Date: Jul 2020

Calibration Code:  512636768133 RSE/CSE Ratio: 12 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 5-0.05 EU/mL Standard Curve Linearity: -0.998

Standard Curve Mean Reaction Times:	<u>5.0</u> EU/mL	<u>126</u> seconds
	<u>0.5</u> EU/mL	<u>343</u> seconds
	<u>0.05</u> EU/mL	<u>767</u> seconds

Archived Spike Concentration: 1.01 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst  Date 29 Aug 2018

Reviewed By  Date 29 Aug 2018

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA