



## Endosafe<sup>®</sup> - PTS<sup>™</sup> Cartridges Certificate of Analysis

Reorder Code: PTS2005F Cartridge Lot #:  6676132 Expiration Date: Jun 2022

Calibration Code:  513839665688 RSE/CSE Ratio: 11 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 5-0.05 EU/mL Standard Curve Linearity: -0.999


Standard Curve Mean Reaction Times: 5.0 EU/mL 138 seconds  
0.5 EU/mL 350 seconds  
0.05 EU/mL 796 seconds


Archived Spike Concentration: 0.76 EU/mL Negative Control: Pass

**This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.**

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

**CAUTION: DO NOT FREEZE THE CARTRIDGES**

Qualified Analyst  Date 26 Jun 2020

Reviewed By  Date 27 June 2020

Charles River Laboratories, Inc.  
1023 Wappoo Road, Suite 43-B  
Charleston, SC 29407 USA