



Endosafe[®] - PTS[™] Cartridges Certificate of Analysis

Reorder Code: PTS20005F Cartridge Lot #:  8847131 Expiration Date: Dec 2020

Calibration Code:  422167645181 RSE/CSE Ratio: 12 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 0.5-0.005 EU/mL Standard Curve Linearity: -0.998

Standard Curve Mean Reaction Times:

<u>0.5</u> EU/mL	<u>221</u> seconds
<u>0.05</u> EU/mL	<u>528</u> seconds
<u>0.005</u> EU/mL	<u>1076</u> seconds

Archived Spike Concentration: 0.071 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst T. Bon Myer Date 18 Dec 2018

Reviewed By [Signature] Date 20 Dec 2018

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA