

## **CERTIFICATE OF ANALYSIS**

Standardization of CSE against RSE Pyrotell® Gel-Clot Method

ACC Internal Control Number: ICQ-0622-119 Test date: 29 JUN 2022

Reagent	Sensitivity (EU/mL)	Catalog No.	Lot No.	Expiration Date (YYYY-MM-DD)
Pyrotell <sup>®</sup>	0.03	G5003	522-06-022	2027-06-27
Control Standard Endotoxin (CSE), 0.5 μg/vial (= 500 ng/vial)		E0005	180	2025-11-19
Reference Standard Endotoxin (RSE)		N/A	H0K354	N/A

Potency of CSE (per ng)	<b>10 EU/ng</b> or 10 IU/ng	
Potency of CSE (per vial)	10 EU/ng x 500 ng/vial = <b>5,000 EU/vial</b> or 5,000 IU/vial	

If the CSE is reconstituted with 5 mL, the endotoxin concentration will be 1,000 EU/mL.

MARL	29JUN 2022
Quality Control	Date
Horle lle A	13 Jul 2022
Quality Assurance	Date
Kewi Catty	n. 13 JUL 2022
Quality Assurance	Date

## Notes:

- 1. The stated potency applies only to the combination of CSE and LAL lots specified on this certificate.
- 2. The equivalence of the EU and IU is stated in the USP Bacterial Endotoxins Test chapter, <85>, and in the EP Bacterial Endotoxins chapter, 2.6.14.
- 3. Potency is subject to change within the error of the test at any time.

M\_COA\_QC\_0085

Rev. 5, December 2021

FERA\_616