

Endosafe® - PTS™ Cartridges Certificate of Analysis

Reorder Code: <u>PTS20005F</u> Cartridge Lot #: 6622164 Expiration Date: May 2022
Calibration Code: 422867346985 RSE/CSE Ratio: 9 EU/ng RSE Lot #: H0K354
Archived Standard Curve Range: <u>0.5-0.005</u> EU/mL Standard Curve Linearity: <u>-0.999</u>
Standard Curve Mean Reaction Times: 0.5 EU/mL 228 seconds 0.05 EU/mL 518 seconds 0.005 EU/mL 1073 seconds
Archived Spike Concentration: 0.089 EU/mL Negative Control: Pass
This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.
Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.
CAUTION: DO NOT FREEZE THE CARTRIDGES
Qualified Analyst / San y Date o 4 M A 4 2020 Reviewed By Annual Date of may 2020
Charles River Laboratories, Inc. 1023 Wappoo Road, Suite 43-B

Charleston, SC 29407 USA