



**Endosafe[®] - PTS[™]/MCS[™] Inhibition/Enhancement Screening Cartridges
Certificate of Analysis**

Reorder Code: PTS220

Expiration Date: Nov 2019



Calibration Code: 512034043539



Cartridge Lot #: 7826808

Archived Spike Concentration: 0.55 EU/mL

Positive Product Control Results: Reactive/Pass

This lot of PTS/MCS Cartridges has been tested and meets Quality Control testing requirements for positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst *[Signature]*

Date 30 NOV 2017

Reviewed By *[Signature]*

Date 30 Nov 2017

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA