

Endosafe® - PTS™ Cartridges Certificate of Analysis

| Reorder Code: PTS2001F Cartridge Lot #: 1773147 Expiration Date: Sep 2023 |
|--|
| Calibration Code: 119757142478 RSE/CSE Ratio: 10 EU/ng RSE Lot #: H0K354 |
| Archived Standard Curve Range: <u>1-0.01</u> EU/mL Standard Curve Linearity: <u>-0.999</u> |
| Standard Curve Mean Reaction Times: 1.0 EU/mL 197 seconds 0.1 EU/mL 469 seconds 0.01 EU/mL 971 seconds |
| Archived Spike Concentration: 0.088 EU/mL Negative Control: Pass |
| This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results. |
| Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only. |
| CAUTION: DO NOT FREEZE THE CARTRIDGES |
| Qualified Analyst Date Date |
| Reviewed By C. Line Date of Nov 2021 |
| Charles River Laboratories, Inc. 1023 Wappoo Road, Suite 43-B Charleston, SC 29407, USA |