

Endosafe[®] - PTS[™] Cartridges
Certificate of Analysis

Reorder Code: PTS2005F Cartridge Lot #:  8718138 Expiration Date: Aug 2020

Calibration Code:  512334978595 RSE/CSE Ratio: 12 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 5-0.05 EU/mL Standard Curve Linearity: -0.999

Standard Curve Mean Reaction Times: 5.0 EU/mL 123 seconds
0.5 EU/mL 325 seconds
0.05 EU/mL 749 seconds

Archived Spike Concentration: 1.05 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst JZJ Date 09 Aug 2018
Reviewed By A. Anderson Date 09 Aug 2018

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA