

CERTIFICATE OF ANALYSIS

VIAL CONTENTS: E contains 10 ng of purified Lip	indosafe [®] Control Standa opolysaccharide, freeze d	rd Endotoxin is pre Iried in a stabilized	epared fi matrix.	rom <i>E. coli</i> strain	055:B5.	Each vial
RSE/CSE RATIO: Th	e potency of this standard or Kinetic Techniq u		, (EU) ha	as been determin	ed to be	
CSE Lot: EX149	052 CSE Exp. Date:	Apr 2024	-			
LAL Reagent Lot:N3343	BL LAL Exp. Date:	May 2025	-			
RSE Lot: H0K3	54 RSE/CSE Ratio:	18	_EU/ng	Vial contents: _	180	_ EU/vial
IS/CSE RATIO: The Expert Committee on Biological Standardization of WHO has assigned a potency of the IS as 10,000 IU per vial of IS, so that 1 IU = 1 EU. The potency of this endotoxin standard in International (Endotoxin) Units, IU, has been designated as18IU/ng. DIRECTIONS FOR USE: Reconstitute the lyophilized material with3.6mL of LAL reagent grad water to obtain50EU/mL or50IU/mL. Vortex mix vigorously for at least 5 minutes after rehydration, and for at least 1 minute immediately prior to each use.						
STORAGE: Store rehydrated material at 2-8°C for up to 4 weeks. Store lyophilized material at controlled room temperature or refrigerated as preferred. Diluted endotoxin should not be stored except under validated conditions.						
CAUTION: DO NOT FREEZE ENDOTOXIN SOLUTIONS						
Prepared By:	Mil		Dat	e: <u>05 Jan</u> e: <u>05 Jan</u>	2022	-
Quality Reviewed By	: Brooke Mart	me	Dat	e: 05 Jan	२०२२	3