



Endosafe[®] - PTS[™] Cartridges Certificate of Analysis

Reorder Code: PTS201F Cartridge Lot #:  8568128 Expiration Date: Mar 2020

Calibration Code:  612936496685 RSE/CSE Ratio: 9 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 10-0.1 EU/mL Standard Curve Linearity: -0.998

Standard Curve Mean Reaction Times:

<u>10.0</u> EU/mL	<u>129</u> seconds
<u>1.0</u> EU/mL	<u>349</u> seconds
<u>0.1</u> EU/mL	<u>764</u> seconds

Archived Spike Concentration: 1.72 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst  Date 20 Mar 2018

Reviewed By  Date 22 Mar 2018

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA