

Endosafe® - PTS™ Cartridges Certificate of Analysis

Reorder Code: PTS20005F Cartridge Lot #: 1638141 E	xpiration Date: <u>May 2023</u>
Calibration Code: 422467046742 RSE/CSE Ratio: 11 EU/ng	RSE Lot #: <u>H0K354</u>
Archived Standard Curve Range: <u>0.5-0.005</u> EU/mL Standard Curve	e Linearity:
Standard Curve Mean Reaction Times: 0.5 EU/mL 224	seconds
0.05 EU/mL 532	seconds
0.005EU/mL1070	seconds
Archived Spike Concentration: <u>0.087</u> EU/mL Negative Contr	ol: <u>Pass</u>
This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.	
Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.	
CAUTION: DO NOT FREEZE THE CARTRIDGES	
Qualified Analyst / Son my Date o	0950N2021
Reviewed By Date	15 Jun 2021
	Charles River Laboratories, Inc.

Charleston, SC 29407 USA