



**Endosafe[®] - PTS[™] Cartridges
Certificate of Analysis**



Reorder Code: PTS201F Cartridge Lot #: 7658130 Expiration Date: Jun 2019



Calibration Code: 611532097197 RSE/CSE Ratio: 15 EU/ng RSE Lot #: HOK354

Archived Standard Curve Range: 10-0.1 EU/mL Standard Curve Linearity: -0.997

Standard Curve Mean Reaction Times: 0.1 EU/mL 115 seconds

1.0 EU/mL 325 seconds

10.0 EU/mL 720 seconds

Archived Spike Concentration: 1.82 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst J. Ly

Date 21 Jun 2017

Reviewed By MBain

Date 28 Jun 2017

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA