

## Endosafe® - PTS™ Cartridges Certificate of Analysis

4500447 = 1 4 5 4 4 5 0000
Reorder Code: PTS20005F Cartridge Lot #: 1528147 Expiration Date: Jan 2023
Calibration Code: 421165744198 RSE/CSE Ratio: 9 EU/ng RSE Lot #: H0K354
Archived Standard Curve Range: <u>0.5-0.005</u> EU/mL Standard Curve Linearity: <u>-0.999</u>
Standard Curve Mean Reaction Times:0.5EU/mL211seconds
0.05 EU/mL 508 seconds
0.005 EU/mL 1057 seconds
Archived Spike Concentration: 0.061 EU/mL Negative Control: Pass
This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.
Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened.  Cartridges are for single-test use only.
CAUTION: DO NOT FREEZE THE CARTRIDGES
Qualified Analyst Jhu Date 29 Jan 2021
Reviewed By Date _01 Feb 707.
Charles River Laboratories, Inc. 1023 Wappoo Road, Suite 43-B

Charleston, SC 29407 USA