



## Endosafe<sup>®</sup> - PTS<sup>™</sup> Cartridges Certificate of Analysis



Reorder Code: PTS2005F Cartridge Lot #: 6836141 Expiration Date: Dec 2022



Calibration Code: 515636076723 RSE/CSE Ratio: 9 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 5-0.05 EU/mL Standard Curve Linearity: -0.999

Standard Curve Mean Reaction Times: 5.0 EU/mL 156 seconds

0.5 EU/mL 366 seconds


0.05 EU/mL 760 seconds

Archived Spike Concentration: 0.87 EU/mL Negative Control: Pass


**This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.**

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

**CAUTION: DO NOT FREEZE THE CARTRIDGES**

Qualified Analyst 

Date 04 Dec 2020

Reviewed By 

Date 05 Dec 2020

Charles River Laboratories, Inc.  
1023 Wappoo Road, Suite 43-B  
Charleston, SC 29407 USA