


Endosafe® - PTS™/MCS™ Inhibition/Enhancement Screening Cartridges
Certificate of Analysis

Reorder Code: PTS220

Expiration Date: Aug 2020

Calibration Code: 
511035093316

Cartridge Lot #: 
8741822

Archived Spike Concentration: 0.53 EU/mL

Positive Product Control Results: Reactive/Pass

This lot of PTS/MCS Cartridges has been tested and meets Quality Control testing requirements for positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst *Radjanica Jott*

Date *20 Sep 2018*

Reviewed By *J. Farmer*

Date *20 Sep 2018*

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA