

Endosafe[®] - PTS[™] Cartridges
Certificate of Analysis

Reorder Code: PTS2001F Cartridge Lot #:  8642156 Expiration Date: May 2020

Calibration Code:  117351606947 RSE/CSE Ratio: 10 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 1-0.01 EU/mL Standard Curve Linearity: -0.997

Standard Curve Mean Reaction Times:

<u>1.0</u> EU/mL	<u>173</u> seconds
<u>0.1</u> EU/mL	<u>443</u> seconds
<u>0.01</u> EU/mL	<u>916</u> seconds

Archived Spike Concentration: 0.178 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst  Date 14 Jun 2018

Reviewed By  Date 20 Jun 2018

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA