



Endosafe[®] - PTS[™] Cartridges Certificate of Analysis

Reorder Code: PTS2005F Cartridge Lot #:  6742155 Expiration Date: Aug 2022

Calibration Code:  515439776494 RSE/CSE Ratio: 11 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 5-0.05 EU/mL Standard Curve Linearity: -1.000

Standard Curve Mean Reaction Times:

<u>5.0</u> EU/mL	<u>154</u> seconds
<u>0.5</u> EU/mL	<u>362</u> seconds
<u>0.05</u> EU/mL	<u>797</u> seconds

Archived Spike Concentration: 0.84 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst *IS on ry* Date 24 Sep 2020

Reviewed By *Burman* Date 24 Sep 2020

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA