



ASSOCIATES OF
CAPE COD
INCORPORATED

Specialists in Endotoxin and Glucan Detection

CERTIFICATE OF ANALYSIS
Standardization of CSE against RSE
Pyrotell®-T, Kinetic Turbidimetric Method

ACC Internal Control Number: ICQ-0326-122
Test date: 19 MAR 2026

Reagent	Catalog No.	Lot No.	Expiration Date (YYYY-MM-DD)
Pyrotell®-T	T0051	525-11-059-T	2030-12-08
Control Standard Endotoxin (CSE), 0.5 µg/vial (= 500 ng/vial)	E0005	194	2029-05-30
Reference Standard Endotoxin (RSE)	N/A	R172R0	N/A

Potency of CSE (per ng)	13 EU/ng or 13 IU/ng
Potency of CSE (per vial)	13 EU/ng x 500 ng/vial = 6,500 EU/vial or 6,500 IU/vial
Valid range of potency	0.001 – 4.0 EU/mL

If the CSE is reconstituted with 5 mL, the endotoxin concentration will be 1,300 EU/mL.

Further dilution will be required to obtain a CSE concentration of 1,000 EU/mL.

	24 Mar 26
Quality Control	Date
	26 MAR 26
Quality Assurance	Date
	30 MAR 26
Quality Assurance	Date

Notes:

1. The stated potency applies only to the combination of CSE and LAL lots specified on this certificate.
2. The equivalence of the EU and IU is stated in the USP Bacterial Endotoxins Test chapter, <85>, and in the EP Bacterial Endotoxins chapter, 2.6.14.
3. Potency is subject to change within the error of the test at any time.

M_COA_QC_0681

Rev. 6, June 2023
FERA_775