

Endosafe® - PTS™ Cartridges Certificate of Analysis

Reorder Code: PTS2001F Cartridge Lot #: 6767129	Expiration Date: Sep 2022
Calibration Code: 118861944524 RSE/CSE Ratio: 11	EU/ng RSE Lot #: H0K354
Archived Standard Curve Range: 1-0.01 EU/mL Standa	ard Curve Linearity:
Standard Curve Mean Reaction Times:1.0 _EU/mL0.1 _EU/mL0.01 _EU/mL	461 seconds
Archived Spike Concentration: 0.130 EU/mL Negative	ve Control: Pass
This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.	
Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.	
CAUTION: DO NOT FREEZE THE CARTRIDGES	
Qualified Analyst Reviewed By Bureun	Date <u>6 6 Nov 2020</u>
	Charles River Laboratories, Inc. 1023 Wappoo Road, Suite 43-B Charleston, SC 29407 USA