



Endosafe[®] - PTS[™] Cartridges Certificate of Analysis



Reorder Code: PTS201F Cartridge Lot #: 6640166 Expiration Date: May 2022



Calibration Code: 614233775179 RSE/CSE Ratio: 8 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 10-0.1 EU/mL Standard Curve Linearity: -0.999

Standard Curve Mean Reaction Times: 10.0 EU/mL 142 seconds

1.0 EU/mL 340 seconds

0.1 EU/mL 737 seconds

Archived Spike Concentration: 1.42 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst MEcay

Date 22 May 2020

Reviewed By B. M. Leung

Date 22 May 2020

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA