

## Endosafe® - PTS™ Cartridges Certificate of Analysis

Reorder Code: PTS20005F Cartridge Lot #: 1523166	Expiration Date: Jan 2023
Calibration Code: 421866543547 RSE/CSE Ratio: 11 E	U/ng RSE Lot #:H0K354
Archived Standard Curve Range: <u>0.5-0.005</u> EU/mL Standard	d Curve Linearity:
Standard Curve Mean Reaction Times: 0.5 EU/mL 0.05 EU/mL 0.005 EU/mL	510 seconds
Archived Spike Concentration: 0.055 EU/mL Negative	Control: Pass
This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.	
Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.	
CAUTION: DO NOT FREEZE THE CARTRIDGES	
Qualified Analyst	Date <u>27 Jan 202</u>
promise the second seco	Date 27 Jan 2021
	Charles River Laboratories, Inc.

1023 Wappoo Road, Suite 43-B Charleston, SC 29407 USA