



## Endosafe<sup>®</sup> - PTS<sup>™</sup> Cartridges Certificate of Analysis



Reorder Code: PTS2001F Cartridge Lot #: 9776127 Expiration Date: Oct 2021



Calibration Code: 115544704583 RSE/CSE Ratio: 9 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 1-0.01 EU/mL Standard Curve Linearity: -0.997

Standard Curve Mean Reaction Times:

<u>1.0</u> EU/mL	<u>155</u> seconds
<u>0.1</u> EU/mL	<u>404</u> seconds
<u>0.01</u> EU/mL	<u>847</u> seconds

Archived Spike Concentration: 0.130 EU/mL Negative Control: Pass

**This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.**

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

**CAUTION: DO NOT FREEZE THE CARTRIDGES**

Qualified Analyst MEcay Date 08 Oct 2019

Reviewed By [Signature] Date 09 Oct 2019

Charles River Laboratories, Inc.  
1023 Wappoo Road, Suite 43-B  
Charleston, SC 29407 USA