

## Endosafe® - PTS™ Cartridges Certificate of Analysis

Reorder Code: PTS5505F Cartridge Lot #: 1879121 Expiration Date: Oct 2023
Calibration Code: 512137353678 RSE/CSE Ratio: 10 EU/ng RSE Lot #: H0K354
Archived Standard Curve Range: <u>5-0.05</u> EU/mL Standard Curve Linearity: <u>-0.995</u>
Standard Curve Mean Reaction Times:5.0 _EU/mL121 _seconds
0.5 _EU/mL357 _seconds
0.05 EU/mL 773 seconds
Archived Spike Concentration: 0.56 EU/mL Negative Control: Pass
This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.
Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened.  Cartridges are for single-test use only.
CAUTION: DO NOT FREEZE THE CARTRIDGES
Qualified Analyst <u>Lumbleyn</u> Date <u>OI Dec 2011</u>
Reviewed By Dateoa bec 202
Charles River Laboratories, Inc. 1023 Wappoo Road, Suite 43-B Charleston, SC 29407 USA