



## Endosafe<sup>®</sup> - PTS<sup>™</sup> Cartridges Certificate of Analysis



Reorder Code: PTS20005F Cartridge Lot #: 7777125 Expiration Date: Oct 2019



Calibration Code: 417150105182 RSE/CSE Ratio: 10 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 0.5-0.005 EU/mL Standard Curve Linearity: -0.993

Standard Curve Mean Reaction Times: 0.5 EU/mL 171 seconds  
0.05 EU/mL 463 seconds  
0.005 EU/mL 901 seconds

Archived Spike Concentration: 0.071 EU/mL Negative Control: Pass

**This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.**

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

**CAUTION: DO NOT FREEZE THE CARTRIDGES**

Qualified Analyst *R. M.*

Date 18 Oct 2017

Reviewed By *M. Allen*

Date 18 Oct 2017

Charles River Laboratories, Inc.  
1023 Wappoo Road, Suite 43-B  
Charleston, SC 29407 USA