

## Endosafe® - PTS™ Cartridges Certificate of Analysis

Reorder Code: PTS2005F Cartridge Lot #: 2505164	4 Expiration Date: Jan 2024
Calibration Code: 515137384456 RSE/CSE Ratio: 9	EU/ng RSE Lot #: H0K354
Archived Standard Curve Range: <u>5-0.05</u> EU/mL Standa	ard Curve Linearity:
Standard Curve Mean Reaction Times:5.0EU/mL	372 seconds
Archived Spike Concentration: 0.64 EU/mL Negative	ve Control: Pass
This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.	
Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.	
CAUTION: DO NOT FREEZE THE CARTRIDGES	
Qualified Analyst	Date 12 Jan 2022
Reviewed By	Date 12 Jan 2022
	Charles River Laboratories, Inc. 1023 Wappoo Road, Suite 43-B Charleston, SC 29407 USA