




Endosafe[®] - PTS[™] Cartridges Certificate of Analysis

Reorder Code: PTS20005F Cartridge Lot #:  7643152 Expiration Date: May 2019

Calibration Code:  420167304789 RSE/CSE Ratio: 13 EU/ng RSE Lot #: H0K354

Archived Standard Curve Range: 0.5-0.005 EU/mL Standard Curve Linearity: -0.995


Standard Curve Mean Reaction Times: 0.5 EU/mL 201 seconds
0.05 EU/mL 534 seconds
0.005 EU/mL 1073 seconds


Archived Spike Concentration: 0.067 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst  Date 10 Oct 2017

Reviewed By  Date 10 Oct 2017

Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43-B
Charleston, SC 29407 USA