## BIOENDO

Xiamen Bioendo Technology Co., Ltd.

PASS

endo Technol

CERTIFICATE OF ANALYSIS				
Catalog Number:	KT5017S	QA	Lot Number:	25020120
Product:	Bioendo KT	Endotoxin Test Kit		
(Kinetic Turb		bidimetric Assay), 800 Tests/Kit	<b>Expiration Date:</b>	05-Feb-2027
Contents		Description	Catalog No.	Lot No.
Lyophilized Amebocyte Lysate		1.7ml/vial	KT17S	25024102
Reconstitution Buffer		3.0ml/vial	KTRB30	25028115
Control Standard Endotoxin		Bacterial Strain: E.coli O111:B4	CSE10V	25025108
CSE Reconstituted is 50EU/ml				
CSE Reconstitution Volume (ml)		2.8		
Test (Method)		Specifications		Results
Linearity		r  ≥0.98		Pass
Enzyme Activity (Reaction Time)				
Blank-0.005EU/ml		$\geq$ 120 seconds		Pass
0.005EU/ml		$\leq 5100$ seconds		Pass
1.0EU/ml		$\leq$ 1500 seconds		Pass
Coefficient of Variation (%)		≤ 10%	20	Pass

Reconstitute CSE with 2.8 ml Water for BET. Mix vigorously for 5 minutes on a vortex mixer, obtaining 50 EU/ml standard stock solution. Reconstituted endotoxin is stable for up to 7 days at 2 to 8° C.

The endotoxin level is 140EU/vial when assay parallel with USP Reference Standard Endotoxin Lot No.: H0K354.

Operator

Lechen Yom Buelin Die

Date 2025.02.20 Date 2025.02.20

**QA** Administrator

ManufacturerXiamen Bioendo Technology Co., Ltd.Address131 Xinjia Road,Haicang Xinyang Industry Zone, Xiamen,ChinaEnglish Website: http://www. www.bioendo.com